BCEIA
2005

自动顶空进样器 *AutoHS Headspace Sampler*


COLINTech

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
Chengdu Colin Analysis Technology Ltd.

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Headspace Samplers eliminate tedious and error-producing steps required in other GC sample preparation techniques, enabling you to extract more information from your samples in less time. The intuitive touch screen graphical user interface provides easy control of the sampler. A patent device with dynamic and static function extended the range of the headspace technology application. Reached an advanced level in the world.

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An intuitive, color touch screen user interface provides instant access to system functionality. Set parameters, recall stored methods and start sampling quickly. The modern screen graphical user interface displays a sampling process. And making user know well it.

中英文可选择 Chinese and English are Available

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
This function will help user to complete the whole carrier gas system test with a vial and get a high reliability.


动态和静态任意选择 Dynamic and Static Mode are Available

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Double holes in needle and gas selecting provide dynamic and static headspace methods. Employed the technique prevent from sample decomposition or oxidization. Employing the Purge mode, can sweep the air in vials with carrier gas. The mode is useful in kinetic studies.

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Volatile compounds in almost any sample matrix can be analyzed simply and quickly by headspace gas chromatography. Headspace sampling provides a method that extract sample of liquid and solid with a solvent-free and minimizes rework by eliminating many of the tedious and error-producing. The system completely rely on volatilization to separate the compounds from the sample matrix, so sample extraction, clean up and preconcentration are not necessary.

File Name	Sample Name	Time [min]	Ethanol		
			Area [μs]	Adjusted Area	Height [s]
YC312007.rst	ethanol	1.72977255	0.98	2.3117	
YC312008.rst	ethanol	1.72980135	0.98	2.3138	
YC312009.rst	ethanol	1.72979352	0.98	2.3095	
YC312010.rst	ethanol	1.72977649	0.98	2.3040	
YC312011.rst	ethanol	1.72989283	0.99	2.2975	
YC312012.rst	ethanol	1.72993786	0.99	2.3032	
YC312013.rst	ethanol	1.72993836	0.99	2.3007	
yc312005-20050315-125559.rst	ethanol	1.72991335	0.99	2.3199	
yc312006-20050315-130411.rst	ethanol	1.72987391	0.99	2.3181	
Averages		1.72985558	0.99	2.3087	
%RSD		1.5e-08	1	0.71	
				0.3389	

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AutoHS provides high sample throughput. Equipped with a 40-sample tray and 15 positions in the heating oven allowing unattended overnight operation for many applications. It can load 13 vials to be purged in the Purge Mode. In Constant Mode, a automated constant heating time is employed with optimized sample overlapping for maximum sample throughput in the minimum amount of time. It can simultaneously thermostat up to 12 samples, ensuring that the next sample is always ready for analysis upon completion of the previous run. This adds up to a significant time savings since there is no need to wait for the next sample to heat and get best reproducibility. A deactivated sampling path, from switch valve to transfer line, is standard on the AutoHS. The clean gas swept the sample path.

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In the Progressive Mode, the vials are thermostatted for an increasing time, representing a multiple of the first thermostating time. It is used in the determination of the necessary equilibration time, where a number of vials containing the same sample are to be thermostatted for increasing periods. Automatically optimized samples.

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The technique requires multiple extractions from the sample vial. This is very useful for developing method and running special sample.

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Provided ports to connect with any gas chromatograph with a standard interface. They includes: ReadyIn, StartIn, StopIn, ReadyOut and RunOut.

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
The microcomputer control system is made in special company. The precision pneumatic elements keep the disturbing away for electronics. Gas elements are from other country.

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Address: Room 212, Kexie Building, No.11,
4th Duan, Renminnan Rd. 610041

Tel: +86-28-85260395

Fax: +86-28-85253094