Magnetic ImmunoChromatographic Test MAR™ Assay Development System

MagnaBioSciences

a Quantum Design company

MAR™ 研发分析系统 Product Overview

***Harnessing the proven power and sensitivity
of magnetic induction technology
to enhance your assay development efforts***

**MICT™ ...is the new force
in quantitative immunoassays!**


MICT™ is powered by MAR™ (Magnetic Assay Reader): proven, patented induction technology that provides quantification of biochemical analytes via the detection of superparamagnetic particles (superPMPs) that are bound to the analyte. The detector measures the local magnetic field expressed by the total mass of label in the sample region. By way of a calibration curve based on the mass of magnetic labels in the immune complexes, this value may be correlated to the concentration of analyte in the specimen.

The result? Highly sensitive, accurate assays!

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

Troponin I, Hepatitis B Surface Antigen 与 Staphylococcus Enterotoxin B 标准曲线

Shown below is data from Troponin I, Hepatitis B and Staphylococcal Aureus Enterotoxin assays, in magnetic particle format, indicating both the linearity of the assays (over >2 logs), as well as sensitivity.


Troponin I Standard Curve

HbsAg Standard Curve

Staphylococcal Enterotoxin B

Lateral Flow Assays Incorporating Superparamagnetic Particles

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Lateral flow rapid test chemistry starter kit:

SuperParaMagnetic particles (carboxylated, 200nm)
Conjugation chemistry
Sample pad
Conjugation pad
Nitrocellulose membrane
Wick pad
Top laminate
MICT™ cassettes

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