



**CHARACTERIZATION OF
SALMONELLA TYPHIMURIUM
MEMBRANE PROTEINS
SEPARATED BY 2D GELS**

A X I M A

APPLICATION NOTE

CHARACTERIZATION OF *SALMONELLA TYPHIMURIUM* MEMBRANE PROTEINS SEPARATED BY 2D GELS

INTRODUCTION

Matrix-assisted laser desorption ionization time-of-flight mass spectrometry (MALDI)⁽¹⁾ is now widely used in the field of proteomics. Key requirements in this area are: high sensitivity, high resolution and mass accuracy and rapid analysis for high throughput. Both precursor and product ion (in the form of PSD) data can be used for the positive identification of proteins. Until now, with conventional reflectron designs⁽²⁾, the acquisition of post-source decay spectra (PSD) has been quite time consuming requiring multiple analyses at varied reflectron voltages and stitching together of the resulting spectra. The incorporation of the curved-field reflectron⁽³⁾ into MALDI TOF instruments significantly increased the speed and sensitivity of the PSD analysis by enabling the acquisition of a fully focused PSD spectrum in a single analysis: seamless PSD (sPSD). All the advantages of the curved-field reflectron have now been integrated into a high resolution/high sensitivity MALDI TOF instrument⁽⁴⁾, the AXIMA-CFR.

Proteomics provides a powerful tool for exploring protein levels at the whole cell level. It can be used to determine the response of bacteria to their environment and as such, provides a fingerprint of the status of the cell at a given point in time. This strategy can be used to examine membrane preparations of bacteria for the presence of surface proteins and transport proteins. The AXIMA-CFR was used to analyze a membrane fraction of *Salmonella typhimurium* proteins which had been resolved by 2D-gel electrophoresis (Figure 1).




Figure 1. 2D Gel of the membrane fraction of *Salmonella typhimurium* proteins.

Spots identified by MS analysis:

1. OmpD (outer membrane protein D);
2. FlhC (phase-I Flagellin);
3. OmpA and EF-Tu (outer membrane protein A precursor and elongation factor-Tu)

METHOD

Protein spots were treated by in-gel digestion using trypsin, the peptides extracted and then desalted prior to MALDI sample preparation. The curved-field reflectron MALDI TOF mass spectrometer (AXIMA-CFR) was used, which comprises: UV nitrogen laser (337nm), near ultra high vacuum system with novel laser beam focusing mechanism for increased ion extraction efficiency and increased sensitivity, an integrated 1 GHz transient recorder, optimized ion optics for high resolution and accuracy, a high resolution ion gate for parent ion selection for sPSD analyses, and a microplate-type target.

All the MALDI analyses, for each sample, were performed on a single 0.5 μ l aliquot using the dried droplet sample loading technique with α -cyano-4-hydroxy cinnamic acid matrix saturated in 1:1 acetonitrile/0.1% TFA.

Database searches were performed using the Mascot™ search engine (Matrix Science).

RESULTS

The performance of the instrument in relation to in-gel digest of 2D gel separated proteins is illustrated with an example from *Salmonella typhimurium* (protein spot 3). The following observations were made:-

- 1 A 40 kDa protein (Coomassie Blue stained) digested with trypsin generated a complex series of signals in the mass range 800 to 3500 Da (Figure 2). Isotopic resolution was obtained across the whole mass range at a single pulsed extraction value. This enabled the monoisotopic peak from each signal to be used in the Mascot™ peptide mass fingerprint search engine and resulted in the primary identity of the protein to be outer membrane protein A (OMP A) from *Salmonella typhimurium* (Table 1).
- 2 Fourteen sPSD spectra were generated for 12 of the signals from a single loading of the desalted sample to the target. Selected sPSD spectra were used to confirm the identity of the protein using the Mascot™ MS/MS ion search engine. The results from these searches (Figures 3-7) indicated the presence of a second protein in the 2D gel spot - Elongation factor-Tu (EF-Tu). A review of the parent mass searches indicated that EF-Tu was also present in the search results though at a lower score. A complete analysis of the parent and sPSD spectra is summarized in Table 2 indicating that there is strong evidence for the presence of both proteins in the same 2D gel spot. Assignment of all the signals from the parent spectrum to the sequences of the two proteins indicated that 60% sequence coverage was obtained for OMP A and 30% for EF-Tu.
- 3 sPSD analysis and database searching of the 1641 Da signal (Figure 5) corresponded to a peptide, LGYPITDDLDFYTR, found in a variant of OMP A that was different from the one identified from the parent mass fingerprint. The published sequence of this variant had only one amino acid difference (F114 \leftrightarrow V114) from the form identified by the peptide mass fingerprint. The equivalent peptide, LGYPITDDLDFYTR, was also observed in the parent spectrum, though not subjected to sPSD, which indicated that both forms of OMP A were present in the 2D gel protein spot.
- 4 sPSD analysis of the 1617 Da signal was useful in identifying one of two OMP A peptides of similar mass. sPSD fragmentation confirmed that the signal corresponded to OMP A 250-263 (theoretical mass $MH^+ = 1616.74$) and not OMP A 46-59 (theoretical mass $MH^+ = 1617.79$).
- 5 Further, sPSD analysis of the 1233 Da peak was useful in identifying one of two OMP A peptides of similar mass. sPSD fragmentation confirmed that the signal corresponded to OMP A 326-334 (theoretical mass $MH^+ = 1233.62$) and not EF-Tu 118-128 (theoretical mass $MH^+ = 1233.64$).
- 6 Even low intensity signals such as 1803 Da generated useful sPSD data (Figure 6) for the structural confirmation and correct protein identity.
- 7 Some signals (e.g. 1438 Da) did not correspond to peptides from any of the three proteins identified. These signals may represent post transitional modifications and will be subject to further investigation by sPSD.




Figure 2. Parent spectrum indicating signals from Outer Membrane Protein A and Elongation Factor Tu




Figure 3. sPSD of 1264.4 Da precursor supporting the identity of Outer Membrane Protein A




Figure 4. sPSD of 2304.8 Da precursor supporting the identity of Outer Membrane Protein A




Figure 5. sPSD of 1641.8 Da precursor supporting the identity of Outer Membrane Protein A variant

Accession	Name	Score	Description
Q14333	OUTER MEMBRANE PROTEIN A	44	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	47	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	48	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	49	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	50	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	51	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	52	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	53	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	54	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	55	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	56	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	57	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	58	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	59	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	60	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	61	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	62	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	63	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	64	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	65	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	66	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	67	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	68	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	69	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	70	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	71	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	72	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	73	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	74	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	75	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	76	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	77	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	78	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	79	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	80	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	81	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	82	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	83	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	84	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	85	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	86	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	87	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	88	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	89	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	90	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	91	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	92	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	93	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	94	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	95	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	96	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	97	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	98	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	99	OUTER MEMBRANE PROTEIN A
Q14333	OUTER MEMBRANE PROTEIN A	100	OUTER MEMBRANE PROTEIN A

Table 1. Mascot™ database search results




Figure 6. sPSD of 1803.8 Da precursor supporting the identity of Elongation Factor-Tu




Figure 7. sPSD of 1966.2 Da precursor supporting the identity of Elongation Factor-Tu

CONCLUSIONS

The integration of the curved-field reflectron into a high performance MALDI TOF mass spectrometer has enabled:

- Very high sensitivity: sPSD analysis : ≤ 50 fmol
- Good sequence coverage for database searching
- High speed analysis: 60 to 90 sec per sPSD analysis
- Typical 8 to 10 seamless PSD spectra can be obtained from a single 0.5 μ l loading of a digested Coomassie Blue stained protein and 3 to 4 from a silver stained protein
- The data obtained from sPSD provided critical support for the identification and structural characterization of proteins from 2D gels which frequently contained more than one entity

	OMP A variant 1	OMP A variant 2	EF-Tu
Parent spectrum number of peptides database position ¹	10 1	9 2	8 14
sPSD ² - parent peak Da			
1027			B
1233			B
1264	A		
1381		B	
1438	?	?	?
1471	A		
1617			B
1641		A	
1903			A
1965			A
2303		B	
2616		B	

¹protein identified with Mascot™ peptide mass fingerprinting search engine

²A = protein identified with Mascot™ MS/MS ion search engine

B = protein identity confirmed from the theoretical fragment comparison

Table 2. Summary of proteins identified from seamless PSD spectra

ACKNOWLEDGEMENTS

We would like thank N. Sienkiewicz², A. D. Cronshaw^{1,2}, F. B. Ward² and M. P. Gallagher² (The Wellmet Facility¹, Institute of Cell and Molecular Biology², Biology Division, University of Edinburgh, Mayfield Road, Edinburgh EH9 3JR) for preparing the *Salmonella typhimurium* proteins by 2D gels electrophoresis and *in situ* trypsin digestion.

REFERENCES

1. Tanaka, K. et al; *Rapid Comm. Mass Spectrom.*, **2**, (1988), 151-3
2. Mamyryn, B.A.; USSR patent no 198,034, (1966)
3. Cotter, R.J. et al; *Rapid Comm. Mass Spectrom.*, **8**, (1994), 781-5
4. Bowdler, A. et al; ASMS Conf. Mass Spectrom. Allied Topics (1999)

AXIMA-CFR

UK/Europe

Wharfside, Trafford Wharf Road, Manchester, M17 1GP, UK

Telephone: +44 (0)161 888 4400 Facsimile: +44 (0)161 888 4402 e-mail: maldi.sales@kratos.co.uk

North America

100 Red Schoolhouse Road, Building A, Chestnut Ridge, New York 10977, USA

Telephone: +1 (845) 426 6700 Facsimile: +1 (845) 426 6192 e-mail: maldi@kratos.com

Japan

Shimadzu Corporation, 380-1 Horiyamashita, Hadano-City, Kanagawa 259-1304, Japan

Telephone: +81 463 88 8660 Facsimile: +81 463 88 8670 e-mail: kawabata@shimadzu.co.jp

Germany

Shimadzu Deutschland GmbH, Albert-Hahn Strasse 6-10, D-47269 Duisburg, Germany

Telephone: +49 (203) 768 70 Facsimile: +49 (203) 711 734 e-mail: mr@shimadzu.de

Australia

Shimadzu (Oceania) PTY Ltd, Unit T, 10-16 South Street, Rydalmere, NSW 2116, Australia

Telephone: +61 (2) 684 4200 Facsimile: +61 (2) 684 4055 e-mail: john.hewetson@shimadzu.com.au

Pacific Rim

Shimadzu (Asia Pacific) PTE Ltd, 16 Science Park Drive#01-01, The Pasteur, Singapore Science Park, Singapore 0511

Telephone: +65 778 6280 Facsimile: +65 779 2935 e-mail: regie@shimadzu.com.sg

Website

http://www.kratos.com



Certificate number ...1907/99...