

AXIMA

Characterization of *Salmonella typhimurium* membrane proteins separated by 2D gels

The AXIMA-CFR high performance MALDI-TOF with curved field reflectron, was used to analyze a membrane fraction of *S. typhimurium* proteins resolved by 2D gel electrophoresis.

The rapid acquisition of high quality seamless PSD data allowed unequivocal identification and characterization of co-migrating proteins which upon tryptic digestion had generated peptides of a similar nominal mass.

Characterization of *Salmonella typhimurium* Membrane Proteins Separated by 2D Gels

Introduction

Matrix-assisted laser desorption ionization time-of-flight mass spectrometry (MALDI)⁽¹⁾ is now widely used in the field of proteomics. Key requirements in this area are: high sensitivity, high resolution and mass accuracy and rapid analysis for high throughput. Both precursor and product ion (in the form of PSD) data can be used for the positive identification of proteins. Until now, with conventional reflectron designs⁽²⁾, the acquisition of post-source decay spectra (PSD) has been quite time consuming requiring multiple analyses at varied reflectron voltages and stitching together of the resulting spectra. The incorporation of the curved-field reflectron⁽³⁾ into MALDI TOF instruments significantly increased the speed and sensitivity of the PSD analysis by enabling the acquisition of a fully focused PSD spectrum in a single analysis: seamless PSD (sPSD). All the advantages of the curved-field reflectron have now been integrated into a high resolution/high sensitivity MALDI TOF instrument⁽⁴⁾, the AXIMA-CFR.

Proteomics provides a powerful tool for exploring protein levels at the whole cell level. It can be used to determine the response of bacteria to their environment and as such, provides a fingerprint of the status of the cell at a given point in time. This strategy can be used to examine membrane preparations of bacteria for the presence of surface proteins and transport proteins. The AXIMA-CFR was used to analyze a membrane fraction of *Salmonella typhimurium* proteins which had been resolved by 2D gel electrophoresis (Figure 1).

Method

Protein spots were treated by in-gel digestion using trypsin, the peptides extracted and then desalted prior to MALDI sample preparation. The curved-field reflectron MALDI TOF mass spectrometer (AXIMA-CFR) was used, which comprises: UV nitrogen laser (337 nm), near ultra high vacuum system with novel laser beam focusing mechanism for increased ion extraction efficiency and increased sensitivity, an integrated 1 GHz transient recorder, optimized ion optics for high resolution and accuracy, a high resolution ion gate for parent ion selection for sPSD analyses, and a microplate-type target.

All the MALDI analyses, for each sample, were performed on a single 0.5 µl aliquot using the dried droplet sample loading technique with α-cyano-4-hydroxy cinnamic acid matrix saturated in 1:1 acetonitrile/0.1% TFA.

Database searches were performed using the Mascot™ search engine (Matrix Science).




Figure 1. 2D Gel of the membrane fraction of *Salmonella typhimurium* proteins. Spots identified by MS analysis: 1. OMP D (outer membrane protein D); 2. FliC (phase-I flagellin); 3. OMP A and EF-Tu (Outer Membrane Protein A precursor and Elongation Factor-Tu)

Results

The performance of the instrument in relation to in-gel digest of 2D gel separated proteins is illustrated with an example from *Salmonella typhimurium* (protein spot 3). The following observations were made:-

- 1 A 40 kDa protein (Coomassie Blue stained) digested with trypsin generated a complex series of signals in the mass range 800 to 3500 Da (Figure 2). Isotopic resolution was obtained across the whole mass range at a single pulsed extraction value. This enabled the monoisotopic peak from each signal to be used in the Mascot™ peptide mass fingerprint search engine and resulted in the primary identity of the protein to be outer membrane protein A (OMP A) from *Salmonella typhimurium* (Table 1).

Accession	Mass	Score	Description
1. G1122349	37566	66	OUTER MEMBRANE PROTEIN A PROTEUSIS (OUTER MEMBRANE PROTEIN
2. G1122349	37719	57	OUTER MEMBRANE PROTEIN A PROTEUSIS (OUTER MEMBRANE PROTEIN
3. G1122349	44130	52	probable integrase - <i>Shigella flexneri</i> phage phi-4-57
4. G1122349	26126	49	OUTER MEMBRANE PROTEIN A (OUTER MEMBRANE PROTEIN II)
5. G1122349	26114	40	outer membrane protein A - <i>Escherichia fergusonii</i> ATCC 354
6. G1122349	37170	47	OUTER MEMBRANE PROTEIN A PROTEUSIS (OUTER MEMBRANE PROTEIN
7. G1122349	37168	44	(AF234287) outer membrane protein A (<i>Escherichia coli</i>)
8. G1122349	60737	49	(AB00074) ATP-dependent Clp protease subunit <i>Hyphomicrobium</i>
9. G1122349	42150	43	flagellin - <i>Stenotrophomonas</i>
10. G1122349	42150	42	Chain A, Elongation Factor Complex EF-Tu/EF-Ts From <i>Escherichia</i>
11. G1122349	37073	41	SDSA PROTEIN
12. G1122349	43155	41	Chain A, Intact Elongation Factor From <i>E. coli</i>
13. G1122349	43256	41	ELONGATION FACTOR TU (EF-TU)
14. G1122349	43256	41	ELONGATION FACTOR TU (EF-TU) (P-47)
15. G1122349	43224	41	translation elongation factor EF-Tu.A - <i>Salmonella typhimurium</i>
16. G1122349	43256	41	translation elongation factor EF-Tu.A - <i>Escherichia coli</i>
17. G1122349	107087	41	MINOR TRYPTOPHAN RIGIDITY-DEFICIENT PROTEIN OMP
18. G1122349	37348	41	(A1139078) possible sugar nucleotide epimerase/dehydratase
19. G1122349	12940	41	GTP-BINDING PROTEIN COXA
20. G1122349	107371	40	(AF136191) Beilin (<i>Dactylopusia</i> sp. HB-606)

Accession	Mass	Score	Description				
1. G1122349	37566	67	OUTER MEMBRANE PROTEIN A PROTEUSIS (OUTER MEMBRANE PROTEIN III) (O				
Observed	Mr (exact)	Mr (scale)	Delta	Start	End	Mass	Peptide
872.00	871.07	871.51	-0.44	223	229	0	GVKIKIK
918.20	918.19	918.52	-0.33	98	103	0	AGQVSTAK
1069.10	1062.10	1062.54	-0.44	223	231	0	IVLFLHFK
1287.27	1286.26	1286.89	-0.63	322	333	0	KALIKGLAPD
1293.26	1292.26	1292.43	-0.20	118	128	1	LOGHWKADTE
1264.06	1263.06	1263.68	-0.62	216	227	0	DCYVFLVDTF
1261.15	1260.15	1260.65	-0.50	268	280	0	IGDQATWGLGSRK
1470.93	1469.92	1469.70	0.22	01	34	1	HYFQKQKMGATK
1608.10	1607.10	1607.81	-0.71	104	117	0	LVPTTDLDFPTE
2302.70	2301.29	2302.20	-0.91	195	217	0	YCGQKAFVYAFAPAFAPFQTE

Accession	Mass	Score	Description				
2. G1122349	37719	57	OUTER MEMBRANE PROTEIN A PROTEUSIS (OUTER MEMBRANE PROTEIN II)				
Observed	Mr (exact)	Mr (scale)	Delta	Start	End	Mass	Peptide
872.00	871.07	871.51	-0.44	334	340	1	IVKIKIK
918.20	918.19	918.52	-0.33	98	103	0	AGQVSTAK
1073.27	1076.26	1076.88	-0.62	294	303	1	GIYAFKIKAD
1069.10	1062.10	1062.54	-0.44	224	232	0	IVLFLHFK
1287.27	1286.26	1286.89	-0.63	313	333	0	KALIKGLAPD
1293.26	1292.26	1292.43	-0.18	118	128	1	LOGHWKADTE
1464.09	1460.90	1469.81	0.87	104	117	0	LVPTTDLDFPTE
1803.41	1802.40	1802.98	-0.58	293	299	1	AGQVSTLIGESIKRHK
2672.02	2671.02	2670.29	0.73	257	281	1	IGQVFLVDTFEGDQATWGLGSRK

Accession	Mass	Score	Description				
13. G1122349	43256	42	ELONGATION FACTOR TU (EF-TU)				
Observed	Mr (exact)	Mr (scale)	Delta	Start	End	Mass	Peptide
827.25	826.24	826.49	-0.25	118	124	0	KKILIK
1027.27	1026.26	1026.50	-0.24	271	280	0	AGQVSTLIG
1214.29	1213.28	1213.62	-0.34	308	314	0	FEEVFLIK
1218.34	1217.33	1217.84	-0.51	178	188	0	ALGKDAKRAK
1293.26	1292.26	1292.41	-0.15	324	334	0	GTSPQVTF
1617.19	1616.18	1615.70	0.49	250	263	1	ETVSTCTGVKRF
1803.41	1802.40	1802.88	-0.48	60	78	0	GIYVSTVETVTF
1964.10	1963.17	1962.95	-0.15	156	172	0	ELIKVDFVDFQVTFVTF

Table 1. Mascot™ database search results

2 Fourteen sPSD spectra were generated for twelve of the signals from a single loading of the desalted sample to the target. Selected sPSD spectra were used to confirm the identity of the protein using the Mascot™ MS/MS ion search engine. The results from these searches (Figures 3-7) indicated the presence of a second protein in the 2D gel spot - Elongation Factor-Tu (EF-Tu). A review of the parent mass searches indicated that EF-Tu was also present in the search results though at a lower score. A complete analysis of the parent and sPSD spectra is summarized in Table 2 indicating that there is strong evidence for the presence of both proteins in the same 2D gel spot. Assignment of all the signals from the parent spectrum to the sequences of the two proteins indicated that 60% sequence coverage was obtained for OMP A and 30% for EF-Tu.




Figure 2. Parent spectrum indicating signals from Outer Membrane Protein A and Elongation Factor-Tu

3 sPSD analysis and database searching of the 1641 Da signal (Figure 5) corresponded to a peptide, LGYPITDDLVDYTR, found in a variant of OMP A that was different from the one identified from the parent mass fingerprint. The published sequence of this variant had only one amino acid difference (F114 <-> V114) from the form identified by the peptide mass fingerprint. The equivalent peptide, LGYPITDDLDFYTR, was also observed in the parent spectrum, though not subjected to sPSD, which indicated that both forms of OMP A were present in the 2D gel protein spot.




Figure 3. sPSD of 1264.6 Da precursor supporting the identity of Outer Membrane Protein A

4 sPSD analysis of the 1617 Da signal was useful in identifying one of two OMP A peptides of similar mass. sPSD fragmentation confirmed that the signal corresponded to OMP A 250-263 (theoretical mass $MH^+ = 1616.74$) and not OMP A 46-59 (theoretical mass $MH^+ = 1617.79$).




Figure 4. sPSD of 2304.8 Da precursor supporting the identity of Outer Membrane Protein A

6 Even low intensity signals such as 1803 Da generated useful sPSD data (Figure 6) for the structural confirmation and correct protein identity.

7 Some signals (e.g. 1438 Da) did not correspond to peptides from any of the three proteins identified. These signals may represent post transitional modifications and will be subject to further investigation by sPSD.




Figure 5. sPSD of 1641.8 Da precursor supporting the identity of Outer Membrane Protein A variant

Figure 6. sPSD of 1803.8 Da precursor supporting the identity of Elongation Factor-Tu





Figure 7. sPSD of 1966.2 Da precursor supporting the identity of Elongation Factor-Tu



	OMP A variant 1	OMP A variant 2	EF-Tu
Parent spectrum			
number of peptides	10	9	8
database position ¹	1	2	14
sPSD² - parent peak Da			
1027			B
1233			B
1264	A		
1381		B	
1438	?	?	?
1471	A		
1617			B
1641		A	
1903			A
1965			A
2303		B	
2616		B	

¹protein identified with Mascot™ peptide mass fingerprinting search engine
²A = protein identified with Mascot™ MS/MS ion search engine
B = protein identity confirmed from the theoretical fragment comparison

Table 2. Summary of proteins identified from seamless PSD spectra

Conclusions

The integration of the curved-field reflectron into a high performance MALDI TOF mass spectrometer has enabled:

- Very high sensitivity: sPSD analysis : ≤ 50 fmol
- Good sequence coverage for database searching
- High speed analysis: 60 to 90 sec per sPSD analysis
- Typical 8 to 10 seamless PSD spectra can be obtained from a single 0.5 μ l loading of a digested Coomassie Blue stained protein and 3 to 4 from a silver stained protein
- The data obtained from sPSD provided critical support for the identification and structural characterization of proteins from 2D gels which frequently contained more than one entity

Acknowledgements

We would like thank N. Sienkiewicz², A. D. Cronshaw^{1,2}, F. B. Ward² and M. P. Gallagher² (The Welmet Facility¹, Institute of Cell and Molecular Biology², Biology Division, University of Edinburgh, Mayfield Road, Edinburgh EH9 3JR) for preparing the *Salmonella typhimurium* proteins by 2D gel electrophoresis and in situ trypsin digestion.

REFERENCES

1. Tanaka, K. et al; *Rapid Comm. Mass Spectrom.*, 2, (1988), 151-3
2. Mamyrin, B.A.; *USSR patent no 198,034*, (1966)
3. Cotter, R.J. et al; *Rapid Comm. Mass Spectrom.*, 8, (1994), 781-5
4. Bowdler, A. et al; *ASMS Conf. Mass Spectrom. Allied Topics* (1999)