

Transilluminators UVstar

Filter type and intensity setting	UV filter size (cm x cm)			
	15 x 15	20 x 20	23 x 30	20 x 20 + 20 x 20 White light
High-grade filter				
50 or 100 % dual switch	UVstar 15	UVstar 20 UVstar 20HM UVstar 20ML	UVstar 30 UVstar 30HM UVstar 30ML	UVstar WL
10 – 100 % 10 % increment regulator	UVstar 15i	UVstar 20i	UVstar 30i	UVstar WLi
Super Brilliant high-grade filter				
50 or 100 % dual switch	-	UVstar 20 plus	UVstar 30 plus	-
10 – 100 % 10 % increment regulator	-	UVstar 20i plus	UVstar 30i plus	-

Models in bold red are also available without UV protection shield and are compatible with BDA Box 1/2/3.

Note:

All models are available with 254 nm, 312 nm or 312 nm bulbs, but Super Brilliant models only with 312 nm or 365 nm.

White light transilluminators WLstar

Intensity setting	White light filter size (cm x cm)	
	20 x 20	23 x 30
50 or 100 % dual switch	WLstar 20	WLstar 30