

**Multi Dispenser EDR-500A Determination Results (Fluorescent Dye, FITC-I)**

Model: 96 Channels 1-20ul Cylinder

Volume: 20:5.0ul (7 times average)

Spots	A	B	C	D	E	F	G	H	Average	CV
1	984.00	1065.00	1075.00	1033.00	1071.00	1012.00	1115.00	1092.00	<b>1055.88</b>	<b>3.83%</b>
2	1026.00	1019.00	1097.00	1104.00	1066.00	1077.00	1059.00	1065.00	<b>1064.13</b>	<b>2.65%</b>
3	1045.00	1069.00	1110.00	1079.00	1069.00	1006.00	1036.00	1030.00	<b>1055.50</b>	<b>2.90%</b>
4	1037.00	1053.00	1081.00	1050.00	1001.00	1011.00	1054.00	1057.00	<b>1043.00</b>	<b>2.33%</b>
5	1025.00	1105.00	1090.00	1042.00	1018.00	989.00	996.00	1017.00	<b>1035.25</b>	<b>3.79%</b>
6	968.00	1002.00	1002.00	961.00	1030.00	1023.00	1042.00	985.00	<b>1001.63</b>	<b>2.72%</b>
7	1012.00	1061.00	1036.00	1021.00	1043.00	1012.00	1024.00	983.00	<b>1024.00</b>	<b>2.14%</b>
8	1029.00	1042.00	1072.00	1011.00	1001.00	979.00	960.00	1008.00	<b>1012.75</b>	<b>3.26%</b>
9	991.00	1014.00	1086.00	1097.00	1026.00	1007.00	956.00	991.00	<b>1021.00</b>	<b>4.42%</b>
10	1009.00	1067.00	1064.00	1068.00	1076.00	1072.00	990.00	1037.00	<b>1047.88</b>	<b>2.90%</b>
11	966.00	1029.00	1093.00	1095.00	1060.00	1104.00	1009.00	1019.00	<b>1046.88</b>	<b>4.40%</b>
12	1001.00	1052.00	1090.00	1126.00	1118.00	1044.00	1110.00	1120.00	<b>1082.63</b>	<b>3.92%</b>
<b>Average</b>	<b>1007.75</b>	<b>1048.17</b>	<b>1074.67</b>	<b>1057.25</b>	<b>1048.25</b>	<b>1028.00</b>	<b>1029.25</b>	<b>1033.67</b>	<b>1040.88</b>	<b>1.74%</b>
<b>CV</b>	<b>2.60%</b>	<b>2.74%</b>	<b>2.75%</b>	<b>4.43%</b>	<b>3.31%</b>	<b>3.71%</b>	<b>4.97%</b>	<b>4.16%</b>	<b>2.12%</b>	<b>4.00%</b>

Max.:

1126.00

Min.:

956.00

SD:

41.60